

## Coating procedures for ibidi $\mu$ -Slides and $\mu$ -Dishes

For optimized cell adhesion there are different treatments and coatings for the ibidi labware family. The ibiTreat surface is comparable to standard tissue culture treated labware. This surface permits direct cell growth with a large number of cell lines and primary cells. Compared to the ibiTreat surface, the uncoated surface is very hydrophobic. The uncoated surface must be coated with a protein solution for the adhesion of most cells.

### 1. Recommended Surfaces

For Collagen I: ibiTreat (tissue culture treated)

For Collagen IV: ibiTreat (tissue culture treated) or hydrophobic, uncoated

For Fibronectin: ibiTreat (tissue culture treated) or hydrophobic, uncoated

For Poly-L-Lysine: ibiTreat (tissue culture treated)

For Poly-D-Lysine: ibiTreat (tissue culture treated)

To establish a specific coating relevant to a specific research application, we recommend testing the coating procedure on both uncoated and ibiTreat  $\mu$ -Slides. We have observed that some biomolecules adhere differently to hydrophobic and hydrophilic plastic surfaces.

Please note that there is no uncoated version of the  $\mu$ -Slide Chemotaxis. Use the ibiTreat surface for all coatings.

Some products are also offered with a glass bottom. Glass can also be coated according to the specifications below.

The ESS surface needs higher protein concentrations for an effective coating and must not be dried after coating.

### 2. Prepare the Coating Solution

All coating solutions are calculated for a specific **amount of protein per area** ( $\mu\text{g}/\text{cm}^2$ ) recommended by the manufacturer's reference.

For Collagen I: ( $5 \mu\text{g}/\text{cm}^2$ )

Dilute the Collagen I solution (e.g. ibidi, rat tail, 50202) to the desired concentration using 17.5 mM acetic acid (~0.1% acetic acid).

For Collagen IV: ( $1.5 \mu\text{g}/\text{cm}^2$ )

Dilute the Collagen IV (e.g. Corning, mouse tumor, No. 356233) to the desired concentration using 0.05 M HCl.

For Fibronectin: ( $1.5 \mu\text{g}/\text{cm}^2$ )

Dilute the Fibronectin (e.g. Corning, human plasma, 354008) to the desired concentration using PBS (pH 7.2) without  $\text{Ca}^{2+}$  and  $\text{Mg}^{2+}$ .

For Poly-L-Lysine: ( $2 \mu\text{g}/\text{cm}^2$ )

Dilute the PLL (e.g. Sigma-Aldrich. 0.01% solution, 100  $\mu\text{g}/\text{ml}$ , P4832) to the desired concentration using ultra-pure water.

For Poly-D-Lysine: ( $5 \mu\text{g}/\text{cm}^2$ )

Dilute the PDL (e.g. Corning, No. 354210) to the desired concentration using ultra-pure water.

## Application Note 08

Use the following protein concentrations [ $\mu\text{g/ml}$ ]:

### Channel Slides

|  | Collagen I | Collagen IV | Fibronectin | Poly-L-Lysine | Poly-D-Lysine |
|--|------------|-------------|-------------|---------------|---------------|
| $\mu$ -Slide I                           | 250        | 75          | 75          | 100           | 250           |
| $\mu$ -Slide I 0.2 Luer                  | 500        | 150         | 150         | 200           | 500           |
| $\mu$ -Slide I 0.4 Luer                  | 250        | 75          | 75          | 100           | 250           |
| $\mu$ -Slide I 0.6 Luer                  | 200        | 60          | 60          | 80            | 200           |
| $\mu$ -Slide I 0.8 Luer                  | 125        | 38          | 38          | 50            | 125           |
| $\mu$ -Slide III 3in1                    | 250        | 75          | 75          | 100           | 250           |
| $\mu$ -Slide VI 0.4                      | 250        | 75          | 75          | 100           | 250           |
| $\mu$ -Slide VI 0.1                      | 1000       | 300         | 300         | 400           | 1000          |
| $\mu$ -Slide VI - Flat                   | 250        | 75          | 75          | 100           | 250           |
| $\mu$ -Slide y-shaped                    | 250        | 75          | 75          | 100           | 250           |
| $\mu$ -Slide Chemotaxis <sup>1)</sup>    | 130        | 40          | 40          | 55            | 130           |
| $\mu$ -Slide Chemotaxis <sup>2)</sup>    | 230        | 70          | 70          | 90            | 230           |
| $\mu$ -Slide Chemotaxis 2D <sup>1)</sup> | 150        | 45          | 45          | 60            | 150           |
| $\mu$ -Slide Chemotaxis 2D <sup>2)</sup> | 330        | 100         | 100         | 133           | 330           |
| $\mu$ -Slide Membrane ibiPore Flow       | 250        | 75          | 75          | 100           | 250           |
| $\mu$ -Slide III 3D Perfusion            | 100        | 30          | 30          | 40            | 100           |
| $\mu$ -Slide CorrSight™ Live             | 100        | 30          | 30          | 40            | 100           |

### Open Formats

|   | Collagen I | Collagen IV | Fibronectin | Poly-L-Lysine | Poly-D-Lysine |
|---|------------|-------------|-------------|---------------|---------------|
| $\mu$ -Dish 35 mm, low                    | 50         | 15          | 15          | 20            | 50            |
| $\mu$ -Dish 35 mm, high                   | 50         | 15          | 15          | 20            | 50            |
| $\mu$ -Dish 35 mm, high ESS <sup>3)</sup> | 100        | 30          | 30          | 40            | 100           |
| $\mu$ -Dish 50 mm, low                    | 60         | 18          | 18          | 25            | 60            |
| $\mu$ -Slide 2 well                       | 25         | 8           | 8           | 10            | 25            |
| $\mu$ -Slide 4 well                       | 30         | 9           | 9           | 12            | 30            |
| $\mu$ -Slide 8 well                       | 35         | 11          | 11          | 15            | 35            |
| $\mu$ -Slide 2 well Ph+                   | 38         | 11          | 11          | 15            | 38            |
| $\mu$ -Slide 4 well Ph+                   | 42         | 12          | 12          | 17            | 42            |
| $\mu$ -Slide 2 Well Co-Culture            | 40         | 12          | 12          | 17            | 40            |
| $\mu$ -Slide 18 well - Flat               | 40         | 12          | 12          | 17            | 40            |
| $\mu$ -Slide Angiogenesis                 | 125        | 38          | 38          | 50            | 125           |
| $\mu$ -Plate 24 well                      | 20         | 6           | 6           | 9             | 20            |
| $\mu$ -Plate 96 well                      | 35         | 12          | 12          | 15            | 35            |
| $\mu$ -Plate Angiogenesis 96 well         | 125        | 38          | 38          | 50            | 125           |
| 3 Well Chamber, removable                 | 15         | 5           | 5           | 6             | 15            |
| 8 Well Chamber, removable                 | 35         | 11          | 11          | 15            | 35            |
| 12 Well Chamber, removable                | 35         | 11          | 11          | 15            | 35            |
| Culture-Insert 2 Well                     | 60         | 18          | 18          | 25            | 60            |
| Culture-Insert 3 Well                     | 60         | 18          | 18          | 25            | 60            |
| Culture-Insert 4 Well                     | 60         | 18          | 18          | 25            | 60            |
| micro-Insert 4 Well                       | 115        | 35          | 35          | 47            | 115           |
| micro-Insert 4 Well FulTrac               | 100        | 30          | 30          | 40            | 100           |

<sup>1)</sup> When coating the full chamber.

<sup>2)</sup> When coating the observation area only.

<sup>3)</sup> For the very hydrophobic ESS surface, a higher protein concentration is necessary.

Keep in mind that all walls inside the channels are coated in the ibidi channel slides. Open formats are coated not only on the growth area but also partially on the side walls. The coating areas are valid for the exact coating volumes in the table only.

The dilutions above were calculated using the following coating areas and volumes. The coating area is the area which is in contact with the liquid, thus coated.

### Channel Slides

|                                     | Growth Area [cm <sup>2</sup> ] | Coating Area [cm <sup>2</sup> ] | Coating Volume [µl] |
|-------------------------------------|--------------------------------|---------------------------------|---------------------|
| µ-Slide I                           | 2.5                            | 5.4                             | 100                 |
| µ-Slide I 0.2 Luer                  | 2.5                            | 5.2                             | 50                  |
| µ-Slide I 0.4 Luer                  | 2.5                            | 5.4                             | 100                 |
| µ-Slide I 0.6 Luer                  | 2.5                            | 5.6                             | 150                 |
| µ-Slide I 0.8 Luer                  | 2.5                            | 5.8                             | 200                 |
| µ-Slide III 3in1                    | 1.23                           | 3.05                            | 60                  |
| µ-Slide VI 0.4                      | 0.6 per channel                | 1.2 per channel                 | 30 per channel      |
| µ-Slide VI 0.1                      | 0.17 per channel               | 0.34 per channel                | 1.7 per channel     |
| µ-Slide VI - Flat                   | 0.6 per channel                | 1.2 per channel                 | 30 per channel      |
| µ-Slide y-shaped                    | 2.8                            | 5.6                             | 110                 |
| µ-Slide Chemotaxis <sup>1)</sup>    | 1.24 per chamber               | 3.5 per chamber                 | 130 per chamber     |
| µ-Slide Chemotaxis <sup>2)</sup>    | 0.06 per chamber               | 0.27 per chamber                | 6 per chamber       |
| µ-Slide Chemotaxis 2D <sup>1)</sup> | 0.96 per chamber               | 2.4 per chamber                 | 80 per chamber      |
| µ-Slide Chemotaxis 2D <sup>2)</sup> | 0.07 per chamber               | 0.39 per chamber                | 6 per chamber       |
| µ-Slide Membrane ibiPore Flow       | 1.25 (lower channel)           | 2.7 (lower channel)             | 50 (lower channel)  |
| µ-Slide III 3D Perfusion            | 0.25 per well                  | 2.4 per channel                 | 130 per channel     |
| µ-Slide CorrSight™ Live             | 0.25 per well                  | 2.4 per channel                 | 130 per channel     |

### Open Formats

|                                      | Growth Area [cm <sup>2</sup> ] | Coating Area [cm <sup>2</sup> ] | Coating Volume [µl] |
|--------------------------------------|--------------------------------|---------------------------------|---------------------|
| µ-Dish 35 mm, low                    | 3.5                            | 4.1                             | 400                 |
| µ-Dish 35 mm, high                   | 3.5                            | 4.1                             | 400                 |
| µ-Dish 35 mm, high ESS <sup>3)</sup> | 3.5                            | 4.1                             | 800                 |
| µ-Dish 50 mm, low                    | 7.0                            | 7.9                             | 700                 |
| µ-Slide 2 well                       | 4.8 per well                   | 7.5 per well                    | 1500 per well       |
| µ-Slide 4 well                       | 2.2 per well                   | 4.1 per well                    | 700 per well        |
| µ-Slide 8 well                       | 1.1 per well                   | 2.2 per well                    | 300 per well        |
| µ-Slide 2 well Ph+                   | 4.8 per well                   | 11.4 per well                   | 1500 per well       |
| µ-Slide 4 well Ph+                   | 2.2 per well                   | 5.9 per well                    | 700 per well        |
| µ-Slide 2 Well Co-Culture            | 0.4 per minor well             | 0.55 per minor well             | 70 per minor well   |
| µ-Slide 18 well - Flat               | 0.2 per well                   | 0.25 per well                   | 30 per well         |
| µ-Slide Angiogenesis                 | 0.12 per well                  | 0.23 per well                   | 10 per inner well   |
| µ-Plate 24 well                      | 1.9 per well                   | 4.3 per well                    | 1000 per well       |
| µ-Plate 96 well                      | 0.55 per well                  | 2.35 per well                   | 300 per well        |
| µ-Plate Angiogenesis 96 well         | 0.12 per well                  | 0.23 per well                   | 10 per inner well   |
| 3 Well Chamber, removable            | 1.66 per well                  | 3.37 per well                   | 1100 per well       |
| 8 Well Chamber, removable            | 0.93 per well                  | 2.63 per well                   | 400 per well        |
| 12 Well Chamber, removable           | 0.56 per well                  | 1.9 per well                    | 250 per well        |
| Culture-Insert 2 Well                | 0.22 per well                  | 0.82 per well                   | 70 per well         |
| Culture-Insert 3 Well                | 0.22 per well                  | 0.82 per well                   | 70 per well         |
| Culture-Insert 4 Well                | 0.35 per well                  | 1.23 per well                   | 110 per well        |
| micro-Insert 4 Well                  | 0.03 per well                  | 0.23 per well                   | 10 per well         |
| micro-Insert 4 Well FulTrac          | 0.0012 per well                | 0.188 per well                  | 10 per well         |

<sup>1)</sup> When coating full chamber.

<sup>2)</sup> When coating observation area only.


<sup>3)</sup> Also valid for glass bottom and Grid-50/Grid-500 versions.

Keep in mind all walls inside the channels are coated in the ibidi channel slides. Open formats are coated not only on the growth area but also partially on the side walls. The coating areas are valid for the exact coating volumes in the table only.

**3. Fill the channel or the well with the coating solution using the coating volume from the tables above.**

Quick dispensing helps fill the channel slides more easily. Work under sterile conditions. Incomplete filling, or large air bubbles, leads to reduced coating. The ibiTreat surface is easier to completely wet with the recommended volumes than the hydrophobic, uncoated surface.

The very small channels (channel height 0.2 mm and smaller) are filled more easily by using a small volume syringe with a male Luer tip as shown on the right.




**4. Incubate at room temperature for 60 minutes.**

**5. Aspirate the channel or the well volume completely.**

**6. Rinse carefully with ultra-pure water or PBS.**

For rinsing we recommend using 5-10 times the volume of the channel or well. When rinsing a channel slide you can easily add solution into one channel end and simultaneously aspirate it on the other side as shown.



Rinsing thoroughly is necessary to remove all unbound proteins. Any remaining unbound protein may inhibit cell attachment.

**7. Wells or channels are ready to use. Optionally, let dry at room temperature. **Attention, some coating proteins might degenerate during drying! Coatings on the ESS surface must not be dried!****

**8. Store under sterile conditions and use as soon as possible.**

**IMPORTANT NOTES:**

Due to the fact that adhesion proteins are biological substances, there can be quality differences between the lots of the manufacturer. Therefore, it is recommended to test every lot number prior to large scale experiments. Prepare and use other coating substrates according to the manufacturer's specifications or reference.