# Instructions





LifeAct<sup>®</sup> is a 17 amino acids long fragment of a protein originating from *Saccharomyces cerevisiae*, which comprises an actin–binding domain. This marker can be used in various eukaryotic cells to stain filamentous actin (F–actin). Used in living cells it is perfectly labeling the highly dynamic F–actin and moreover, does not interfere with cellular processes.

### **Vector Description**

p<sup>CAG</sup>–LifeAct<sup>®</sup>–TagRFP is a mammalian expression vector encoding LifeAct<sup>®</sup>–TagRFP fusion protein. The vector can be used for fluorescent labeling of the actin cytoskeleton in various living cells. TagRFP codon usage is optimized for high expression in mammalian cells, i.e. humanized [Haas et al., 1996]. Actin–binding domain of the yeast protein Abp140 is fused to the TagRFP N-terminus [Riedl et al., 2008]. For more information on the reporter please visit www.evrogen.com. p<sup>CAG</sup>–LifeAct<sup>®</sup>–TagRFP vector can be used as a source of LifeAct<sup>®</sup>–TagRFP hybrid sequence. The vector backbone contains unique restriction sites that permit its excision and further insertion into expression vector of choice (XhoI, NotI).

The vector backbone also contains the cytomegalovirus immediate early enhancer coupled to chicken  $\beta$ -actin promoter (CAG) [Niwa et al. 1991] for protein expression and SV40 polyadenylation signals (SV40 poly A) for proper processing of the 3' end of the reporter mRNA. SV40 early promoter (P<sub>SV40</sub>) provides neomycin resistance gene (Neo<sup>R</sup>) expression to select stably transfected eukaryotic cells using G418. Neo<sup>R</sup> gene is linked with herpes simplex virus (HSV) thymidine kinase (TK) polyadenylation signals.



## Specifications

#### **Location of Features**

P<sub>CAG</sub>: 365-1704 CMV IE enhancer: 1-364 TATA box: 616-621

LifeAct<sup>®</sup>: 1717-1767

TagRFP Startcodon: 1789-1791 Stopcodon: 2500-2502

SV40 early mRNA polyadenylation signal Polyadenylation signals: 2655-2660 & 2684-2689 mRNA3'ends: 2690 & 2702

SV40 early promoter Enhancer (72-bp tandem repeats): 3378-3449 & 3450-3521 21-bp repeats: 3525-3545, 3546-3566 & 3568-3588 Early promoter element: 3601-3607

Neomycin resistance gene (Neo<sup>®</sup>) Neomycin phosphotransferase coding sequences: Startcodon: 3732-3734 Stopcodon: 4524-4526

Herpes simplex virus (HSV) thymidine kinase (TK) polyadenylation signal Polyadenylation signals: 4759-4764 & 4772-4777



### **Expression in Mammalian Cells**

p<sup>CAG</sup>–LifeAct<sup>®</sup>–TagRFP can be transfected into mammalian cells by any known transfection method. CAG promoter provides strong, constitutive expression of the LifeAct<sup>®</sup>–TagRFP fusion in eukaryotic cells. If required, stable transformants can be selected using G418 [Gorman, 1985].

### Propagation in E. coli

Suitable host strains for propagation in *E. coli* include DH5alpha, HB101, XL1–Blue, and other general purpose strains. Plasmid incompatibility group is pMB1/ColE1. The vector confers resistance to kanamycin ( $30 \mu g/ml$ ) to *E. coli* hosts. Copy number in *E. coli* is about 500.

#### References

Gorman, High efficiency gene transfer into mammalian cells. In DNA cloning: A Practical Approach, Vol. II. Ed.

D. M. Glover. (IRL Press, Oxford, U.K.), 1985: 143-90

Haas et al., Codon usage limitation in the expression of HIV–1 envelope glycoprotein. Curr Biol, 1996, 6 (3): 315–324

Niwa et al., Efficient selection for high-expression transfectants with a novel eukaryotic vector. Gene, 1991, 108: 193–200

Riedl et al., LifeAct: a versatile marker to visualize F–actin. Nature Methods, 2008, 5 (7): 605–607

#### Note:

The vector sequence has been compiled using the information from sequence databases, and published literature, together with partial sequences obtained by ibidi. This vector has not been completely sequenced.



# Instructions

# **Ordering Information**

# LifeAct<sup>®</sup> Plasmids

| c <sup>x</sup> | Cat. No. | Description   | Amount |
|----------------|----------|---|--------|
| A A            | 60101    | <b>p<sup>CMV</sup>–LifeAct<sup>®</sup>–TagGFP2</b> : plasmid, ready to use, 500 ng/μl | 20 µg  |
| e le           | 60102    | p <sup>CMV</sup> –LifeAct <sup>®</sup> –TagRFP։ plasmid, ready to use, 500 ng/µl      | 20 µg  |
| Yibidi         | 60106    | <b>p<sup>CAG</sup>–LifeAct<sup>®</sup>–TagGFP2</b> : plasmid, ready to use, 500 ng/μl | 20 µg  |
|                | 60107    | p <sup>CAG</sup> –LifeAct <sup>®</sup> –TagRFP: plasmid, ready to use, 500 ng/µl      | 20 µg  |

### LifeAct<sup>®</sup> Adenoviral Vectors

|       | Cat. No. | Description  |            |           |       |    |      | Amount                     |
|-------|----------|--|------------|-----------|-------|----|------|----------------------------|
| Hed o | 60121    | <b>rAV<sup>CMV</sup>–LifeAct<sup>®</sup>–TagGFP2</b> :<br>1 × 10 <sup>10</sup> IU/ml | adenovira  | l vector, | ready | to | use, | $1 \times 10^9  \text{IU}$ |
|       | 60122    |  | adenoviral | vector,   | ready | to | use, | $1 \times 10^9  \text{IU}$ |

## LifeAct<sup>®</sup> Lentiviral Vectors

| ifed c. | Cat. No. | Description   | Amount                 |
|---------|----------|---|------------------------|
|         | 60141    | <b>rLV<sup>Ubi</sup>–LifeAct<sup>®</sup>–TagGFP2</b> : lentiviral vector, ready to use, 1 × 10 <sup>7</sup> TU/ml | 1 × 10 <sup>6</sup> TU |
|         | 60142    | <b>rLV</b> <sup>Ubi</sup> –LifeAct <sup>®</sup> –TagRFP: lentiviral vector, ready to use, $1 \times 10^7$ TU/ml   | 1 × 10 <sup>6</sup> TU |

## LifeAct<sup>®</sup> Cell Lines

| Thidi | Cat. No. | Description  | Amount       |         |            |
|-------|----------|--|--------------|---------|------------|
|       | 40101    | HT-1080 LifeAct <sup>®</sup> –Tag<br>LifeAct <sup>®</sup> –TagGFP2 | GFP2: HT-108 | 0 cells | expressing |



### **Notice to Purchaser**

#### **Evrogen Fluorescent Proteins TagGFP2 and TagRFP**

This product is for internal non–commercial research use only. No rights are conveyed to modify or clone the gene encoding fluorescent protein contained in this product. The right to use this product specifically excludes the right to validate or screen compounds. For information on commercial licensing, contact Evrogen Licensing Department, email: *license@evrogen.com*.

#### LifeAct<sup>®</sup> Product and its Variants

Use of these products and product applications are covered by european and foreign patents and patents pending.

All names containing the designation <sup>®</sup> are registered with the Office for Harmonization in the Internal Market and the U.S. Patent and Trademark Office.

Not–For–Profit Entities: Orders may be placed in the normal manner by contacting your local representative or ibidi Customer Service at +49 (0)89/520 4617 0. At its discretion, ibidi grants Not–For–Profit Entities a non–exclusive, personal, limited license to use this product for non–commercial life science research use only. Such license specifically excludes the right to sell or otherwise transfer this product, its components or derivatives thereof to third parties. No modifications to the protein coding sequence may be made without express written permission from ibidi. Any other use of this product requires a license from ibidi. For–Profit Entities wishing to use this product are required to obtain a license from ibidi. For license information, please contact a licensing representative by email at *licensing@ibidi.com*.

#### For research use only!

Further technical specifications can be found at www.ibidi.com. For questions and suggestions please contact us by e-mail *info@ibidi.de* or by telephone +49 (0)89/520 4617 0. All products are developed and produced in Germany. © ibidi GmbH, Am Klopferspitz 19, 82152 Martinsried, Germany.