# Instructions






LifeAct<sup>®</sup> is a 17 amino acids long fragment of a protein originating from *Saccharomyces cerevisiae*, which comprises an actin–binding domain. This marker can be used in various eukaryotic cells to stain filamentous actin (F–actin). Used in living cells it is perfectly labeling the highly dynamic F–actin and moreover, does not interfere with cellular processes.

## **Vector Description**

p<sup>CMV</sup>–LifeAct<sup>®</sup>–TagGFP2 is a mammalian expression vector encoding LifeAct<sup>®</sup>–TagGFP2 fusion protein. The vector can be used for fluorescent labeling of the actin cytoskeleton in various living cells. TagGFP2 codon usage is optimized for high expression in mammalian cells, i.e. humanized [Haas et al., 1996]. Actin–binding domain of the yeast protein Abp140 is fused to the TagGFP2 N-terminus [Riedl et al., 2008]. For more information on the reporter please visit www.evrogen.com. p<sup>CMV</sup>–LifeAct<sup>®</sup>–TagGFP2 vector can be used as a source of LifeAct<sup>®</sup>–TagGFP2 hy-

## Specifications



TagG	FP2 Fluorescence
Ex.max	483 nm
Em. <sub>max</sub>	506 nm
Find more informatio	n on www.evrogen.com.
Packa	aging and Storage
Amount	$20\mu g$ dissolved in $40\mu l$ TE
Concentration	500 ng/µl
Shipping conditions	+2 - 8°C
Storage conditions	-20°C*

Under proper storage conditions

as indicated on vial.

\*Avoid repeated freeze and thaw cycles.

brid sequence. The vector backbone contains unique restriction sites that permit its excision and further insertion into expression vector of choice (XhoI, NotI).

The vector backbone also contains immediate early promoter of cytomegalovirus ( $P_{CMVIE}$ ) for protein expression and SV40 polyadenylation signals (SV40 poly A) for proper processing of the 3' end of the reporter mRNA. SV40 early promoter ( $P_{SV40}$ ) provides neomycin resistance gene (Neo<sup>R</sup>) expression to select stably transfected eukaryotic cells using G418. Neo<sup>R</sup> gene is linked with herpes simplex virus (HSV) thymidine kinase (TK) polyadenylation signals.

#### **Location of Features**

P<sub>CMVIE</sub>: 1-589 Enhancer region: 59-465 TATA box: 554-560

LifeAct<sup>®</sup>: 619-669

TagGFP2 Startcodon: 691-693 Stopcodon: 1406-1408

SV40 early mRNA polyadenylation signal Polyadenylation signals: 1560-1565 & 1589-1594 mRNA3'ends: 1598 & 1610

SV40 early promoter Enhancer (72-bp tandem repeats): 2286-2357 & 2358-2429 21-bp repeats: 2433-2453, 2454-2474 & 2476-2496 Early promoter element: 2508-2514

Neomycin resistance gene (Neo<sup>®</sup>) Neomycin phosphotransferase coding sequences: Startcodon: 2637-2639 Stopcodon: 3429-3431

Herpes simplex virus (HSV) thymidine kinase (TK) polyadenylation signal Polyadenylation signals: 3667-3672 & 3680-3685

Shelf life



## **Expression in Mammalian Cells**

p<sup>CMV</sup>–LifeAct<sup>®</sup>–TagGFP2 can be transfected into mammalian cells by any known transfection method. CMV promoter provides strong, constitutive expression of the LifeAct<sup>®</sup>–TagGFP2 fusion in eukaryotic cells. If required, stable transformants can be selected using G418 [Gorman, 1985].

### Propagation in *E. coli*

Suitable host strains for propagation in *E. coli* include DH5alpha, HB101, XL1–Blue, and other general purpose strains. Plasmid incompatibility group is pMB1/ColE1. The vector confers resistance to kanamycin ( $30 \mu g/ml$ ) to *E. coli* hosts. Copy number in *E. coli* is about 500.

#### References

Gorman, High efficiency gene transfer into mammalian cells. In DNA cloning: A Practical Approach, Vol. II. Ed. D. M. Glover. (IRL Press, Oxford, U.K.), 1985: 143–90

Haas et al., Codon usage limitation in the expression of HIV–1 envelope glycoprotein. Curr Biol, 1996, 6 (3): 315–324

Riedl et al., LifeAct: a versatile marker to visualize F–actin. Nature Methods, 2008, 5 (7): 605–607

#### Note:

The vector sequence has been compiled using the information from sequence databases, and published literature, together with partial sequences obtained by ibidi. This vector has not been completely sequenced.



# Instructions

## **Ordering Information**

## LifeAct<sup>®</sup> Plasmids

C <sup>*</sup>	Cat. No.	Description	Amount
A	60101	<b>p<sup>CMV</sup>–LifeAct<sup>®</sup>–TagGFP2</b> : plasmid, ready to use, 500 ng/μl	20 µg
e le	60102	p <sup>CMV</sup> –LifeAct <sup>®</sup> –TagRFP։ plasmid, ready to use, 500 ng/µl	20 µg
Yibidi	60106	$p^{CAG}$ –LifeAct <sup>®</sup> –TagGFP2: plasmid, ready to use, 500 ng/µl	20 µg
	60107	p <sup>CAG</sup> –LifeAct <sup>®</sup> –TagRFP: plasmid, ready to use, 500 ng/µl	20 µg

## LifeAct<sup>®</sup> Adenoviral Vectors

	Cat. No.	Description						Amount
le A	60121	<b>rAV<sup>CMV</sup>–LifeAct<sup>®</sup>–TagGFP2</b> : 1 × 10 <sup>10</sup> IU/ml	adenovira	l vector,	ready	to	use,	$1 \times 10^9  \text{IU}$
Ybidi	60122		adenoviral	vector,	ready	to	use,	$1 \times 10^9  \text{IU}$

## LifeAct<sup>®</sup> Lentiviral Vectors

	Cat. No.	Description	Amount
S A	60141	<b>rLV</b> <sup>Ubi</sup> –LifeAct <sup>®</sup> –TagGFP2: lentiviral vector, ready to use, $1 \times 10^7$ TU/ml	$1 \times 10^{6}  \mathrm{TU}$
Weine	60142	<b>rLV<sup>Ubi</sup>–LifeAct<sup>®</sup>–TagRFP</b> : lentiviral vector, ready to use, $1 \times 10^7$ TU/ml	$1 \times 10^{6} \mathrm{TU}$
Tbld			

## LifeAct<sup>®</sup> Cell Lines

<u>,</u>	Cat. No.	De
Tibidi	40101	HT Lif

Cat. No.	Description				Amount	
40101	<b>HT-1080 LifeAct<sup>®</sup>–TagGFP2</b> : LifeAct <sup>®</sup> –TagGFP2	HT-1080	cells	expressing	$5 \times 10^5$ cell	



### **Notice to Purchaser**

#### **Evrogen Fluorescent Proteins TagGFP2 and TagRFP**

This product is for internal non–commercial research use only. No rights are conveyed to modify or clone the gene encoding fluorescent protein contained in this product. The right to use this product specifically excludes the right to validate or screen compounds. For information on commercial licensing, contact Evrogen Licensing Department, email: *license@evrogen.com*.

#### LifeAct<sup>®</sup> Product and its Variants

Use of these products and product applications are covered by european and foreign patents and patents pending.

All names containing the designation <sup>®</sup> are registered with the Office for Harmonization in the Internal Market and the U.S. Patent and Trademark Office.

Not–For–Profit Entities: Orders may be placed in the normal manner by contacting your local representative or ibidi Customer Service at +49 (0)89/520 4617 0. At its discretion, ibidi grants Not–For–Profit Entities a non–exclusive, personal, limited license to use this product for non–commercial life science research use only. Such license specifically excludes the right to sell or otherwise transfer this product, its components or derivatives thereof to third parties. No modifications to the protein coding sequence may be made without express written permission from ibidi. Any other use of this product requires a license from ibidi. For–Profit Entities wishing to use this product are required to obtain a license from ibidi. For license information, please contact a licensing representative by email at *licensing@ibidi.com*.

#### For research use only!

Further technical specifications can be found at www.ibidi.com. For questions and suggestions please contact us by e-mail *info@ibidi.de* or by telephone +49 (0)89/520 4617 0. All products are developed and produced in Germany. © ibidi GmbH, Am Klopferspitz 19, 82152 Martinsried, Germany.